

Laboratory Bottles with Thread, DURAN®, DURAN Group / IDL

Laboratory bottles with graduation, ISO 4796-1, with Retrace Code (Batch Certificate of Performance via Internet). Clear glass according to USP/EP and JP.

① Without cap and pouring ring

Volume ml	Thread GL	O.D. mm	Height mm	Order No.	PU pc/Pack
10 ²⁾	25	36	50	10 21 10108	10
25 ²⁾	25	36	70	11 21 80114	10
50	32	46	87	11 21 80117	10
100 ¹⁾	45	56	100	11 21 80124	10
150 ¹⁾	45	62	110	10 21 10129	10
250 ¹⁾	45	70	138	11 21 80136	10
500 ¹⁾	45	86	176	11 21 80144	10
750 ¹⁾	45	95	203	10 21 10153	10
1000 ¹⁾	45	101	225	11 21 80154	10
2000 ¹⁾	45	136	260	11 21 80163	10
3500	45	160	295	10 21 10162	1
5000	45	182	330	11 21 80173	1
10000	45	227	410	11 21 80186	1
15000	45	268	445	11 21 80188	1
20000	45	288	505	11 21 80191	1



② With screw cap and pouring ring made of PP (blue)

Volume ml	Thread GL	O.D. mm	Height mm	Order No.	PU pc/Pack
10 ²⁾	25	36	50	10 21 15108	10
25 ²⁾	25	36	74	11 21 85114	10
50	32	46	91	11 21 85117	10
100 ¹⁾	45	56	105	11 21 85124	10
150 ¹⁾	45	62	114	10 21 11129	10
250 ¹⁾	45	70	143	11 21 85136	10
500 ¹⁾	45	86	181	11 21 85144	10
750 ¹⁾	45	95	207	10 21 11153	10
1000 ¹⁾	45	101	230	11 21 85154	10
2000 ¹⁾	45	136	265	11 21 85163	10
3500	45	160	299	10 21 11162	1
5000	45	182	335	11 21 85173	1
10000	45	227	415	11 21 85186	1
15000	45	268	450	11 21 85188	1
20000	45	288	510	11 21 85191	1



Accessories

Silicone Bottle Holder

The new silicone bottle holder helps to stabilize bottles during such laboratory activities such as dispensing liquids with bottle mounted dispensers or connection caps, pipetting or filtration. The centre of mass of the laboratory glass bottle is set centrally in the bottle holder, so you can work comfortably and safely. The bottle is held securely in the holder by eight flexible ribs, making it suitable for use with all round and square bottles that have a diameter/edge length of approx. 75 – 120 mm. In addition, the holder can be used for holding other laboratory glassware, such as tall and low beakers or filtration flasks.

Article	Outer Ø mm	Order No.
④ DURAN® Silicone Bottle Holder	165	10 21 15401

1) A bead is molded in the glass on the shoulder of the 100 to 2000 ml size to indicate the nominal capacity level.

2) With a specially formed glass rim for better pouring (hence no need for an additional pouring ring made of plastic).

Important: Autoclave these flasks only with a loosely screwed cap (do not close tightly).

Membrane cap for all autoclavable processes see page XXX. The bottle can be tightly closed with the membrane cap. The membrane provides for essential pressure compensation.

