

Special-Connector for Rotary Evaporators, Rettberg

Adapter for Rotary Evaporator for distillation direct from the storage jar without the need to transfer.

The unit consists of a central funnel tube NS 29/32 which feeds into the shaft using a double-ended PVDF screw section with an internal FPM O-ring to fit storage jars from 5 to 50 ml.

Connector for rotary evaporators

Article	Order No.
Adapter for 5 and 10 ml threaded jars	13 70 04401
Adapter for 15, 20, 30, 40 and 50 ml- threaded jars	13 70 04402

Screw jars for the adapters

Article	Order No.	PU pc/Pack
Screw-top jar 5 ml	13 70 04407	100
Screw-top jar 10 ml	13 70 04408	100
Screw-top jar 15 ml	13 70 04411	100
Screw-top jar 20 ml	13 70 04412	100
Screw-top jar 30 ml	13 70 04416	100
Screw-top jar 50 ml	13 70 04419	100

