## Autoloop PRO – Fully Automated Sterilisation for Inoculation Loops, WLD-TEC

## Suitable for Fuego Series and Gasprofi micro Series

The carrousel rotates and controls flaming automatically and detects removed and replaced inoculation loop. Continuous working during the flaming and the cooling phases of the inoculation loops saves up to 70% time.

- Fully graphic display & stable housing
- Removal position on both sides for right and left handers
- Adjustable flaming and cooling time exact to the second
  During the flame sterilizing, the display shows the remaining flaming time
- Additional flame sterilisation (e.g. of bottle necks), possible with optional foot pedal
- An intelligent sensor monitors the safe sterilization
- Display protected by heat resistant glass
- 4 inoculation loop holders & pedestal

Laboratory gas burner not included.

## **Technical data:**

left/right, depends on removed holder Rotation direction: 1 – 4, with optical holder detection Inoculation loop holder:

1 – 15 s, adjustable Time of flame sterilizing:

3 setting options / off switch Adjustable sterilization monitor: Safety features: temperature control

flame activation warning indicator of remaining cool-down

tilt sensor and residual heat display Stand by automatic unit switch off: 1 h stainless steel, aluminum, UV and solvent resistant

Casing: Measurements: Height Autoloop Pro 305 mm

	Article	Order No.
1	Autoloop Pro for Fuego series, pedestal 135 x 180 mm	10 30 78038
	Autoloop Pro for Gasprofi micro series, pedestal 135 x 135 mm	10 30 78039
	Splash guard for Autoloop Pro	32 30 69021



