SALE UP TO Auto loop PRO **Speed-up Flame Sterilizing Inoculation Loops** and Cut Down Your Costs! Case 1: Bunsen Burner & 1 Inoculation Loop Loop 1 4 sequences finished after 80 seconds Working with inoculation loop Flaming inoculation loop Time to exchange loops Cooling period Case 2: Fuego SCS & 4 Inoculation Loops placed on a Tray Loop 1 Loop 2 Loop 3 Loop 4 ... after 40 seconds Case 3: Fuego SCS with Automated Caroussel **Autoloop Pro & 4 Inoculation Loops** Loop 1 (Optimized Flaming Loop 2 with Sterilization Monitor) Loop 3 Loop 4 after 24 seconds Working with the Autoloop Pro (Case 3): • Save 70% working time compared to Case 1 • Save 40% working time compared to Case 2 Compared to Case 2 the Autoloop Pro has paid off after approx. 50 working hours!*