

TICKOPUR Cleaning Agents, BANDELIN



BANDELIN
Ultraschall seit 1955

Bio-degradable cleaning and disinfection products especially for use with **SONOREX ultrasonic bathes**, gentle on items being cleaned and the cleaning tank, enhances ultrasonic cavitation, can be thinned with water, high yield. Choice of alkali, neutral or acidic, emulsion or non-emulsifying TICKOPUR cleaning agents depending on the cleaning task. 1-l and 2-l bottles, 5-l and 25-l cans. 200-l barrels and other cleaning and disinfecting products on request.

TICKOPUR R 33 - Universal cleaner, EXAM examined

Concentrate, dosage 3 – 5%, mild alkaline, pH 9.9 (1%), also suitable for metal alloy, gentle on material, provides corrosion protection. May be corrosive on parts made of pewter and zink. Can be used for dip and wipe cleaning and for cleaning laboratory equipment and fixtures. Can be used as contact fluid in an ultrasonic bath.

Volume	Order No.	Volume	Order No.
1 l bottle	35 07 16001	5 l canister	35 07 16002
① 2 l bottle	35 07 16003	25 l canister	35 07 16005

TICKOPUR R 30 - Neutral cleaning agent

Concentrate, dosing 1 – 5%, neutral, pH 7, suitable for light metals. Pre-testing is recommended for sensitive materials. Emulsifying, gentle acting, with corrosion protection.

Volume	Order No.	Volume	Order No.
1 l bottle	35 07 15001	5 l canister	35 07 15002
② 2 l bottle	35 07 15003	25 l canister	35 07 15005

TICKOPUR R 27 - Special intensive cleaning agent

Concentrate, based on phosphoric acid, dosing 5%, acidic, pH 1.9 at 1%. Removes lime deposits, mineral residues and rust. Not for light metals, pewter and zinc. Special care is required in case of defective chromium plating.

Volume	Order No.	Volume	Order No.
1 l bottle	35 07 14001	5 l canister	35 07 14002
③ 2 l bottle	35 07 14003	25 l canister	35 07 14005

TICKOPUR RW 77 - Special cleaning agent with ammonia

Concentrate, dosing 5%, alkaline, pH 9.9 at 1%, emulsifying, with ammonia, phosphate-free. Powerful cleaning action with slightly deoxidizing effect. Special care is required in case of light metals.

Volume	Order No.	Volume	Order No.
1 l bottle	35 07 12001	5 l canister	35 07 12002
④ 2 l bottle	35 07 12003	25 l canister	35 07 12005

STAMMOPUR 24 - Desinfection and intensive cleaning-certified (VAH/EXAM)

Free of aldehydes, chlorine and phenolene. Bactericidal, levurocidal, virucidal against vaccines, BVDV, HBV, HCV, HIV, H5N1. Application 1 – 2%, mildly alkaline, pH 9.4 (1%). With corrosion protection. Very high material compatibility.

Volume	Order No.	Volume	Order No.
⑤ 2 l bottle	35 07 19002	25 l canister	35 07 19025
5 l canister	35 07 19005		

TICKOPUR TR 2 - demulsifying - special cleaning agent

Concentrate, dosing 0.1 – 5%, acidic, pH 3.6 at 1%, with corrosion protection. Also removes mineral contaminants, greases, oil and pigments. Special care is required in case of light metal, pewter and zinc.

Volume	Order No.	Volume	Order No.
⑥ 2 l bottle	35 07 10103	25 l canister	35 07 10105
5 l canister	35 07 10102		

TICKOPUR TR 7 - demulsifying universal cleaning agent

Concentrate, dosing 0.1 – 5%, mildly alkaline, pH 8.9 at 1%, with corrosion protection, phosphate-free. Can also be used in immersion and wiping processes.

Volume	Order No.	Volume	Order No.
⑦ 2 l bottle	35 07 10003	25 l canister	35 07 10005
5 l canister	35 07 10002		

TICKOPUR TR 13 - demulsifying - special cleaning agent

Concentrate, dosing 0.1 – 10%, alkaline, pH 11.9 at 1%, not for alloys, tin, zinc. Silicate-free. Removes gumming, soot, coking pigments.

Volume	Order No.	Volume	Order No.
⑧ 2 l bottle	35 07 10203	25 l canister	35 07 10205
5 l canister	35 07 10202		

TICKOPUR TR 3 - special cleaning agent

Concentrate, based on citric acid. Application 5%, weakly acidic, pH 3 at 1%. For non-ferrous, precious, and light metals. Removes mineral residues, flash rust, greases, oils, and waxes.

